




Innovative Hygiene.



LUNA handDESINFECT dispenser	
Product Description Housing made of antistatic, impact resistant ABS plastic, white polished. With elbow operation, refill indicator, lockable.	
Field of application/use LUNA handDESINFECT dispenser with elbow operation, therefore no contact with contaminated hands. The trouble-free changeable vacuumBAG shrink bottle and the dispense tube build the airsealed tank. The operation of the dispenser and the quality of it's contents remain unchanged in case the dispenser hasn't been used for a longer time. With each elbow operation a portion of LUNA handDISINFECT is dispensed.	
Dimension H x B x D mm	270 x 120 x 140
Wrapping units LUNA handDESINFECT dispenser	Art.-No. 005013
Assembling Either above or next to the wash basin, as the dispenser is drip-free. Recommended user height at approx. 1300 mm from the bottom of the dispenser to the floor. In case of mounting above the washbasin, min. 150mm distance between the bottom of the dispenser and the edge of the basin.	
Refills LUNA handDESINFECT, 6 x 700 ml	Art.-No. 00626
References Hagleitner offers with the LUNA dermaPROTECT System an innovative solution for all kinds of working areas. Consisting of LUNA handPROTECT, LUNA handDESINFECT, LUNA abrasivSEIFEN and LUNA handCREAM Spender.	
Package Unit	12
Europalette	360

LUNA handDESINFECT Spender (posuda)	
Opis proizvoda Kućiče od antistatičke ABS-plastike otporne na udare, sjajnobelo. Rukovanje laktom, ima prozorče za kontrolu sadržaja, može se zaključati.	
upotreba/područje upotrebe LUNA handDESINFECT Spender (posuda) pokreće se pomoću lakta i zbog toga nema kontakta sa zaraženim rukama. Lako zamenjiva smanjujuća flaša (vacuumBAG) zajedno sa prednjim crevom, čini rezervoar koji ne propušta vazduh. Funkcija rezervoara kao i kvalitet sadržaja ostaju isti i prilikom dužeg neupotrebljavanja. Prilikom svakog pritiska laktom izlazi jedna porcija LUNA handDESINFECT sredstva.	
Veličina H x B x T mm (V x Š x D = visina x širina x dubina)	270 x 120 x 140
Količina pakovanja LUNA handDESINFECT Spender (posuda)	Br. -Art. 005013
Montaža Po želji iznad ili pored umivaonika zato što spender (posuda) ne kapa. Preporučena visina ca. 1300 mm od donjeg ruba kućišta spendera (posude) do poda. Prilikom montaže iznad umivaonika minimalni razmak između donjeg ruba spendera (posude) i umivaonika 150 mm.	
Za nadopunjavanje LUNA handDESINFECT, 6 x 700 ml	Br. -Art. 00626
Preporuka, ocena LUNA dermaPROTECT System-om Hagleitner nudi inovativno rešenje za sva područja rada. Sastoji se od LUNA handPROTECT, LUNA handDESINFECT, LUNA abrasivSEIFEN und LUNA handCREAM Spender (posuda).	
Pakovanje	12
Evropska paleta	360